
Sorptions-Isotherm

Accessory for *LabMaster-aw* System

INDEX

ENGLISH

How to use the NOVASINA accessory "Sorptions Isotherm Set"	5
Recording of an SI, Sorptions Isotherm curve	6
Example.....	7

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How to use the NOVASINA accessory "Sorpton Isotherm Set"

Water in Foods

Most materials and all food products contain water: bound water and free water.

The WATER or MOISTURE CONTENT is the **weight** in % of the *total* water in a product. Test methods are eg. the Karl Fischer Titration, or drying balances.


The WATER ACTIVITY, a_w , is the free water, the only value, which informs about the **micro-biological safety** of a product, specially important for foodstuffs.

In materials it is called the RELATIVE EQUILIBRIUM HUMIDITY (% reh).

The range of 0,00...1,00 a_w corresponds to 0...100% reh.

The relation between moisture content and water activity is the SI, Sorpton Isotherm curve, which is *product and temperature specific*. An absolutely constant temperature during all measurements is therefore very important!

Examples of SI-curves:

SI - Sorption Isotherm curves of foods
(also called water vapour isotherms)

Recording of an SI, Sorption Isotherm curve

best with a Novasina aw instrument with well sealing measuring chamber and full temperature control, selectable from 0°C to 50°C \pm 0.2K.


Adsorption curve: in order to start, the substance must be dry, so it first has to be totally dried out in an oven (depending on the material eg. 2 hours at 102°C). Then after weighing the dry substance (net weight) on a precision balance, the first possible salt to use is the SAL-T6. After about one hour or more, note the aw value and weigh the substance again. When waiting longer, the aw value as well as the weight are higher, moving always along the SI-curve, so the "equilibration time" is not critical.

Then continue with SAL-T11, up to SAL-T98.



Desorption curve: the substance must be more humid when starting to weigh (to be safely on the desorption curve to avoid SI-hysteresis mistakes), so the first salt to use is SAL-T98. When reaching a value well above 0.9, then SAL-T90 can be placed. After about one hour or more, note the aw value and weigh the substance. When waiting longer, the aw value as well as the weight are lower, moving always along the SI-curve, so the waiting time is not critical.

After that, place SAL-T75, down to SAL-T6. Then, after totally drying the substance in the oven, the net (dry) weight is known and all points can be plotted.

1. Weigh the tare (empty substance basket and weighing glass). 2. Place the large cup insert.

3. Place the first SAL-T standard 4. Fill basket with substance and place on top of salt

5. Close the cover
 6. Select the desired chamber temperature
 7. Wait until the chamber temperature is reached, then wait for the equilibrium (stability observation time factor set as desired).
-

Example

of doing a Sorption Isotherm curve (adsorption) with a **biscuit**.

- Selected temperature at the Novasina instrument: 25,0°C.
- The totally dried out sample has a net weight of 0.485g.
- SAL-T 11% was placed together with the sample in the chamber. After about one hour, the aw-value was 0.087. Weighing the sample showed 0.519g, which is 6.94% more weight; this is the moisture content at this point!
- Next salt: SAL-T 33%. Displayed values: 0.290aw, 0.530g (+ 9.1%).
- Next salt: SAL-T 53%. Displayed values: 0.523aw, 0.553g (+13.84%).
- Next salt: SAL-T 75%. Displayed values: 0.749aw, 0.605g (+24.55%).

Higher humidity was in this case of a "dry biscuit" not of interest.
With these 5 points, an SI-curve can roughly be drawn.

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Novasina, a division of Axair Ltd.