
full line

PH ELECTRODE

PE-06HD, PE-04HD, PE-03K7, PE-05T









The Art of Measurement

full line

pH ELECTRODE

PH ELECTRODE	* General purpose, laboratory and field usage.
Model: PE-03	* 12.3 mm dia. x 160 mm.
	* Epoxy body, 1 - 13 pH.
PH ELECTRODE	* General purpose, laboratory and field usage.
Model: PE-11	* 10 mm dia. x 130 mm.
	* Epoxy body.
	* 1 - 13 pH. (0 - 14 pH typical)
PH ELECTRODE	* Professional, laboratory and field usage.
	* 9.5 mm dia. x 130 mm.
	* Epoxy body.
	* 0 - 14 pH.
SPEAR TIP	* The "Spear Tip pH electrode " is perfect for those pH
PH ELECTRODE	measurements in applications where sample piercing
	is required. Meat, sausage and cheese are ideal
Model: PE-06HD	applications. The electrodefeatures a very durable
PE-04HD	glass measuring spear packaged in a rugged virtually
	unbreakable epoxy body.
	* Range : 0 to 14 pH (PE-04HD)
	* Range : 1 to 13 pH (PE-06HD)
pH ELECTRODE + Temp. PROBE	* One kit combine of PH electrode PE-03 and Temp.
(combination kit of pH electrode)	probe TP-07 as a whole unit.
	* Easy operation for the pH measurement under
Model: PE-03K7	the ATC (Automatic Temp. Compensation).
	* Available for PH-208, PH-207, PH-207HA, PH-221
	YK-2005WA.
	pH electrode build in Temp. thermister inside.
pH electrode build in Temp.	* pH meter when make the temperature compensation,
thermister inside.	* the extra optional Temp. probe is not necessary.
	Available for PH-207, PH-208, PH-207HA, PH-221,
Model: PE-05T	YK-2005WA,
	General purpose, laboratory & field usage. 1 to 13 pH
	* (typical 0 to 14 pH).
	pH Operation Temp. : 5 to 60 $^{\circ}\mathrm{C}$ (41 to 140 $^{\circ}\mathrm{F}$).
	* Temp. Probe Range : 0 to 60 $^{\circ}\mathrm{C}$ (32 to 140 $^{\circ}\mathrm{F}$).
	* Epoxy Body Material.
	* Connector :
	* BNC plug for pH, ear phone plug for Temp.

^{*} Appearance and specifications listed in this brochure are subject to change without notice. 0609-PHELECTRODE